



Modern microwave digestion techniques

Tomi Kemppinen



MILESTONE
HELPING
CHEMISTS

WHAT CAUSE MOST ERRORS IN CHEMICAL ANALYSIS?




Fig. 1 Causes of error in chemical analysis.

Anal. Methods, 2013, 5, 2914

TRACE METALS ANALYSIS INSTRUMENTATION


Trace metals analysis can be performed by:

- AAS; GF-AA
- ICP; ICP-MS

There is a significant transition to ICP-MS (trace analysis) due to:

- New regulation (e.g. USP)
- Demand of lower detection limits
- Higher productivity (than GFAA)

THE BOTTLENECK



- Modern analytical equipment gives an accurate analysis in few minutes/seconds. The conventional sample preparation requires: time, reagent consumptions and handling
- Analytical measurement is as good as the sample preparation

WHY MICROWAVE CLOSED-VESSEL DIGESTION


- Speed/ productivity
 - Microwave closed vessel digestion allow fast digestion of multiple samples simultaneously
- Digestion quality
 - Microwave closed-vessel allow to work at high T&P enhancing the digestion quality
- Safety
 - Microwave hardware ensure high level of safety
- Ease of use
 - Built-in application libraries, easy handling
- Full control
 - Temperature and Pressure control

THE KEY PARAMETERS IN SAMPLE PREPARATION

- Sample type
- Sample weight
- Acid chemistry
- Digestion conditions
(temperature and pressure)



SAMPLE TYPE




SAMPLE AMOUNT vs PRESSURE


- The sample amount directly impact the maximum pressure
- Vessels and rotors must be able to handle the pressure generated
- Releasing pressure during the run could lead to:
 - Loss of volatile elements
 - High amount of acid fumes in the cavity
 - Releasing pressure is not in compliance with official methods

ORGANIC SAMPLE REACTIVITY

Medium reactivity
(low mass of food samples < 0,5 g)



High reactivity
(high mass of food samples > 4 g)



THE EFFECT OF THE TEMPERATURE IN A DIGESTION PROCESS

Test conditions

Test 1

Test 2

Sample

Baby food

Baby food

Sample amount

0,5g

0,5g

Acid mixture

10mL of HNO₃ 65%

10mL of HNO₃ 65%

Vessel Volume

100mL

100mL

Digestion time

20 minutes

20 minutes


Set Temperature

180°C
32 bar

200°C
45 bar

THE EFFECT OF THE TEMPERATURE IN A DIGESTION PROCESS

1. Incomplete digestion due to low T and P (180°C and 32 bar)
2. Complete digestion thanks to the high T and P (200°C and 45 bar)



DIGESTION GUIDELINES

	Typical samples	Typical acid mixture	Typical working T	Comments
Environmental	Soils, sediments, waste water	HNO ₃ or HNO ₃ and HCl	175°C-180°C	Generates low pressure
Food, beverage and feed	Fruits, meat, juices, cereals, butter	HNO ₃ , H ₂ O ₂ , HCl	180°C-200°C	Often reactive samples, generate high P, HCl used for some elements
Specialty chemicals	Polymers, additives	HNO ₃ , H ₂ O ₂ , HCl	200°C-250°C	Challenging samples, reactive and often requires high T to achieve complete digestion
Pharmaceutical	API's, drugs	HNO ₃ , H ₂ O ₂ , HCl	180°C-230°C	Often reactive samples, generate high P, HCl used for some elements
Petrochemical/ Energy	Coal, coke, crude oil, heavy oil	HNO ₃ , H ₂ O ₂ , HCl	200°C-250°C	Challenging samples, reactive and often requires high T to achieve complete digestion
Geological	Ores, rocks	HNO ₃ , HCl, H ₃ PO ₄ , HF	200°C-260°C	Low P, complex sample composition requires various acids
Metals, alloys	Precious metals, catalysts	Aqua regia, HNO ₃ , HCl, HF	200°C-280°C	Often requires high T for long time

TARGET OF SAMPLE PREPARATION

1. Achieve complete digestion and recovery of the target elements based on your analytical equipment (application & performance)
2. Ensure to match your productivity needs (throughput)
3. Easy implementation of the system into your lab workflow (ease of use)
4. Ensure low running cost (ROI & reliability)
5. Get premium support (Expertise and know-how)





QUALITY ANALYSIS
STARTS WITH
GREAT SAMPLE PREP

ETHOS UP

High Performance
Microwave Digestion System

TEMPERATURE CONTROL

- Temperature is the most important parameter in a digestion process
- Temperature affects the digestion quality
- The temperature control enables to simplify the digestion process
- Full temperature control ensures complete digestion and higher safety




EXOTHERMAL REACTIONS

- Often exothermal reactions is noticed during the digestion of organic samples, such as:
 - Food, feed & beverage
 - Pharmaceutical
 - Biological
 - Petrochemical
- System has to be able to control the reaction to ensure:
 - High digestion quality
 - High safety of the process

MILESTONE easyTEMP


- Direct contactless sensor temperature control
- Temperature control in all vessels
- Quality control of the digestion process
- Complete data traceability
- Enhanced safety



easyTEMP CONTROL

- The combination of Milestone easyTEMP sensor and easyCONTROL software ensures full controls of the digestion process
- Integrated PID control
- Power is driven by the temperature control

Digestion of 0,5 g infant formula with 5 mL of HNO₃ 65%



PERFORMANCE AND PRODUCTIVITY

- The digestion rotors directly impact the performance and the productivity in a digestion process
- The choice of the rotor has to take in account:
 - T & P required to achieve complete digestion
 - Productivity
 - Ease of use
 - Reliability of the components



SK-15
High Pressure
Rotor



MAXI-24 HP
High Performance &
Throughput Rotor



MAXI-44
High Throughput Rotor

VENT-AND-RESEAL TECHNOLOGY

MILESTONE PATENT

- Reactive samples and high masses of organic sample, generate an increment of the internal pressure of vessels
- Vent-and-reseal is a safety mechanism to gently release the pressure only in case of overpressure
- This technology improves safety conditions by eliminating the need for a pressure sensor
- Vent-and-reseal ensures:
 - Full recovery of all elements
 - No loss of volatiles
 - High safety



MILESTONE SK-15

- High Pressure and Temperature rotor (260°C, 100 bar)
- 100 mL high purity PTFE vessels
- Suitable for :
 - Hard-to-digest sample
 - High sample amount
 - Highly reactive samples (exothermal reaction)



SK-15 TYPICAL APPLICATION



Agriculture



Beverage



Ceramics



Chemicals



Clinical



Cosmetics



Environmental



Food & Feed



Geochemistry



Metals



Petrochemical



Officials



Pharmaceutical



Pigments



Polymers

MILESTONE MAXI-24 HP

HIGH PERFORMANCE AND THROUGHPUT ROTOR


- High Performance and throughput (260°C, 60 bar)
- 80 mL high purity PTFE vessels
- Suitable for:
 - Moderate to high reactivity samples
 - Moderate to high sample amount
 - High throughput on several matrices



MAXI-24 HP TYPICAL APPLICATION



Agriculture



Beverage



Environmental



Cannabis



Clinical



Food & Feed



Officials



Pharmaceutical

MILESTONE MAXI-44 HIGH THROUGHPUT ROTOR

- High Throughput rotor (260°C, 35 bar)
- 100 mL high purity PTFE vessels
- Suitable for:
 - Easy and medium-to-digest samples
 - Low reactivity (environmental samples)
 - Moderate sample amount



MAXI-44 TYPICAL APPLICATION



Agriculture



Beverage



Environmental



Food & Feed



Officials



Pharmaceutical

EXPERTISE AND KNOW-HOW



easyCONTROL SOFTWARE

- Icons-driven
- Multi-language
- 21CFR- Part 11 compliance
- Built-in application libraries
- More than 300 Application notes
- PDF creator




METHOD DEVELOPMENT

- Milestone library
 - Microwave programs
 - Sample weight and Acid mixture
- Milestone Connect
 - Application notes
 - Application reports
 - Scientific paper library
- Our experience at your service
 - Contact your local Milestone specialist




FLEXIBILITY



INTRODUCING THE SINGLE REACTION CHAMBER (SRC)

- SRC is a newer approach to microwave closed-vessel digestion
- Overcomes most limitations of rotor-based systems




HOW DOES IT WORK?




PRINCIPLE OF OPERATIONS

- Simple vials containing samples and acid are loaded into the chamber





PRINCIPLE OF OPERATIONS

- Nitrogen is automatically loaded into the chamber




PRINCIPLE OF OPERATIONS

- Microwave program begins
- Sample digestion is completed



PRINCIPLE OF OPERATIONS

- Chamber is automatically cooled and vented



PRINCIPLE OF OPERATIONS

- Solutions are ready for dilution and analysis




TEMPERATURE AND PRESSURE IN THE SRC TECHNOLOGY

- SRC enables superior digestion conditions:
 - Up to 300°C
 - Up to 199 bar
- It leads to:
 - Higher digestion quality
 - Complete recovery (no venting, sample charring, or loss of volatiles)
 - Reduced interferences
 - Larger sample masses
 - Complete digestion of challenging samples



EFFICIENT HEATING AND COOLING

- Shorter digestion runs
 - Faster heating
 - Faster cooling by water chiller (12-14 min from 280°C to 80°C)
- Chamber always at low temperature
 - No acid diffusion
 - No overheating
- World's first system with water-cooled magnetron
 - Superior heating efficiency
 - Longer lifetime
 - Noiseless operation



ULTRAWAVE 3 CONSTRUCTION

- High-pressure lines made of:
 - Corrosion-resistant steel
- Separated high-pressure lines:
 - Inlet for gas pre-pressurization
 - Outlet for venting at the end of the run
 - Lower blanks
 - Longer lifetime of the lines



HIGHEST PRODUCTIVITY

- Number of samples per run
- Multiple sample types simultaneously
- Cleaning time



RACKS AND VIALS


- Inexpensive construction
- Vials available in:
 - PTFE-TFM
 - Quartz
 - Disposable Glass

# of position	Volume (mL)	Vials material
7	40	PTFE-TFM, Quartz, Disposable glass
20	15	
27	8	
40	4.5	Glass




SAME TEMPERATURE IN ALL VIALS

- Rack with vials sits in a water base load
 - Homogenous heating
 - Homogenous temperature
 - Optimum heat dissipation in case of exothermal reaction




SAME PRESSURE THROUGHOUT THE ENTIRE CHAMBER

- Pre-pressurization with nitrogen
 - Ensures all the analytes stay in solution
 - Avoids reaching the boiling points of the acids
 - No charring
 - Pressure equilibrium between the chamber and the vials



H₂O BOILING POINT TEMPERATURE VS. PRESSURE



ALL SAMPLE AT ONCE



Any sample, any chemistry, any volume: IN A SINGLE RUN!

Questions!